

TRUFA: User's manual

June 26, 2015

Contents

1	Getting an account	2
2	Usage limitations	2
3	What is this web server for ?	3
3.1	Quick background on RNA-seq analysis	3
3.2	What can be done with TRUFA ?	4
4	Input:	5
5	Output:	6
5.1	Overview:	6
5.2	Statistics provided:	6
5.2.1	Cleaning step:	7
5.2.2	Assembly/Mapping:	7
5.2.3	Identification/Annotation:	8
5.2.4	Expression	8
6	Other useful programs:	9
7	Quick Start:	9
7.1	Launch a Job	9
7.1.1	Upload demo files	10
7.1.2	Start the job	11
7.2	Check output	12
7.2.1	Check reads quality:	12

7.2.2	Check the assembly:	12
7.2.3	Check reads mapping against transcripts	13
7.2.4	Check the Blast2GO results:	13
7.2.5	Check the eXpress results:	13
8	Running an analysis:	14
8.1	Type of input:	14
8.2	Cleaning step:	15
8.3	Assembly and Mapping step:	15
8.4	Identification step:	16
8.5	Expression quantification step:	17
9	Security and privacy note:	17

1 Getting an account

To request an account, send an email to res_support@ifca.unican.es. A link to this email is as well provided on the homepage of [TRUFA](#). Please provide your personal information (first and last name, web page if available) and information about your institution (name, web page). An email with your account information will be sent to you shortly.

If you find TRUFA useful, please cite:

Kornobis, Cabellos, Aguilar, Frias-Lopez, Rozas, Marco & Zardoya (2015). TRUFA: A user-friendly web server for de novo RNA-seq analysis using cluster computing. *Evolutionary Bioinformatics*. 11:97-104.

[Link to article](#)

And acknowledge IFCA at the University of Cantabria for the use of Altamira supercomputer.

2 Usage limitations

In order to provide the best experience to the community of users, these features will be implemented in TRUFA's next release:

- **Each job folder will be kept for 1 month after finalization of the job.** Afterwards, they will be automatically deleted. The process

is currently semi-automatic and will be fully automatized in TRUFA next version.

- **Please limit yourself to 3 simultaneous jobs.** This will be enforced in the next releases fo TRUFA.
- **The Blast against nr and Blast2GO steps are the most computationnaly demanding.** Please perform those steps only on your final assemblies.

3 What is this web server for ?

“RNA-seq, also called whole-transcriptome shotgun sequencing, refers to the use of high-throughput sequencing technologies for characterizing the RNA content and composition of a given sample” ([Wolf, 2013](#)). Until now, analyzing RNA-seq data remains a Bioinformatic challenge. This web server is designed to make analysis of RNA-seq data in a fast and user-friendly manner, by using cluster computing and reducing the amount of Bioinformatic knowledge necessary.

3.1 Quick background on RNA-seq analysis

Recently and especially with the development of new massively parallel sequencing method, the cost of sequencing whole transcriptomes fell drastically. With a thousand dollars, you can get Gigas worth of transcriptome sequences.

In the case of Illumina sequencing, you will receive as output 1 (single-end) or 2 (paired-end) files per sample in fastq format (.fastq, .fq or compressed fastq such as fq.tar.gz, fq.gz ...).

In order to get high throughput, the original RNA sequences have been shredded before sequencing and so your fastq files are filled with millions to billions of sequences, between 50 and 250 bp long each depending on the Illumina sequencing technology. These sequences are called “reads”.

The first step in the RNA-seq analysis will be to clean the reads. Various cleanings can be performed, such as removing duplicates, adapters, poor quality bases/reads, putative contaminants etc...

The second step of the analysis is to reassemble the shredded pieces, the reads, into contigs (named transcripts in the case of RNA-seq) in order to reconstruct the original RNAs. To do so, two main methods are available:

- With a reference genome available, you will use “mapping” methods.
- Without a reference genome: you will use “*de novo* assembly” methods.

So far, TRUFA can only perform *de novo* assembly using the Trinity software. But more is to come ...

Once you obtained the assembly, the next step is to identify those contigs. This is particularly interesting if you are looking for specific genes. For example, if you are interested by venom toxins you could download a venom protein sequences database and then blast your assembly against it. You'll get hits, representing mRNA which could represent precursors of venom proteins. Another way to identify sequences is to scan your assembly for protein profiles (using HMMER for example).

Recently, a specific vocabulary called “Gene Ontology” has been developed in order to help in the identification process. The idea is to link biologically meaningful keywords to nucleotide sequences to identify their biological role, molecular function and location in cellular components. This can be performed using Blast2GO.

Finally, the last main step of an RNA-seq analysis will be to quantify the expression. One of the interest of having high throughput data is that the number of reads sequenced per contig should reflect the expression levels of the corresponding gene. To quantify the expression, the first step will be to align all reads back to the transcripts (using BOWTIE for example). Then, programs such as eXpress can provide you with expression values (such as TPM, FPKM/RPKM or effective counts).



For more basics on RNA-seq

[Wolf \(2013\)](#)

[Martin and Wang \(2011\)](#)

3.2 What can be done with TRUFA ?

So far, TRUFA is allowing you to perform the following steps (programs used are specified in the parentheses):

- Reads cleaning:
 - Quality control (FastQC)
 - Quality trimming and duplicates removal (Prinseq)
 - Trimming adapters (Cutadapt)
 - Filtering out potential contaminants (Blat)
- De novo assembly of your reads (Trinity)
- Reads mapping (Bowtie and Bowtie2)
- Contigs (i.e transcripts) identification:
 - based on sequence alignment (Blat, Blast)
 - based on protein profiles (HMMER)
 - Annotation with GO terms (Blast against nr and Blast2GO)
- Expression quantification:
 - providing TPM, RPKM/FPKM and counts (eXpress, RSEM)

All the steps of the pipeline can be run as a whole or independently, depending on your needs. In the case of making a complete pipeline, your reads files will be cleaned, then assembled into transcripts and finally the transcripts will be identified and quantified.

4 Input:

The input files can be of multiple sources, depending on the kind of analysis you want to perform and the data you have available. So far the accepted inputs are:

- 1 or 2 Illumina reads files (i.e. single-end or paired-end reads) in fastq or compressed (extensions .tar.gz or .gz) fastq format
- 1 file with already assembled reads, i.e contigs in fasta format

**Important!**

On top of the necessary inputs, you will need to specify some important information related to the library construction:

- Size of the insert (for Trinity and Bowtie)
- Adapters sequences (for Cutadapt)
- Strand specific reads ? (for Trinity and Bowtie)

5 Output:

5.1 Overview:

So far TRUFA is providing the classical outputs of a RNA-seq analysis. For a complete analysis, you will get (type of file precised between parenthesis):

- Prinseq report of reads quality (html)
- FastQC report of reads quality (zip and html)
- Reads files after complete cleaning (fastq)
- Assembly file (fasta)
- File with alignment of the reads against the transcripts (bam)
- Results of the Blast2GO annotation (dat)
- Transcript expression quantification (txt)

5.2 Statistics provided:

Following is the list of statistics and graphs obtained in TRUFA depending on the programs used (location in the file manager is indicated in parenthesis):

5.2.1 Cleaning step:

- FASTQC: Total number of reads, quality and GC content per base and per sequence, per base N content, sequence duplication levels, list of over-represented sequences, Kmer content.
(in jobs/Job_X/STAT/fastqc_report/*.fq_fastqc.zip)
- PRINSEQ:
 - Total number of reads and bases, sequence length distribution, GC content distribution, base quality distribution, Occurrence of N, Poly-A/T tails, tag sequence check, sequence duplication, sequence complexity, dinucleotide odds ratios.
(in jobs/Job_X/STAT/prinseq_report/*.stat.html).
 - Trimming report.
(in jobs/Job_X/STAT/*_trim.fq.log.txt).
 - Duplication report.
(in jobs/Job_X/STAT/*_dup.fq.log.txt)

5.2.2 Assembly/Mapping:

- TRUFA's scripts:
 - Transcript length distribution
(in jobs/Job_X/STAT/assembly_qc/Length_distribution.png).
 - Total number of contigs, Length distribution, Total number of bases in the assembly, N50, contigs in N50, GC content.
(in jobs/Job_X/STAT/assembly_qc/assembly_stats.txt).
- BLAST and Trinity script: BlastX hit report against Swiss-Prot, bins for histogram of hits coverage.
(in jobs/Job_X/STAT/assembly_qc/blast_qc/*).
For more information see:
http://trinityrnaseq.sourceforge.net/analysis/full_length_transcript_analysis.html
- CEGMA: Percentage of highly conserved core eukaryotic genes completely and partially covered in the assembly. This gives an estimation of the completeness of the transcriptome sequenced and assembled.
(in jobs/Job_X/STAT/assembly_qc/cegma/*.completeness_report)

- BOWTIE2: Alignment report with percentage of concordant/discordant paired reads and unpaired reads mapped, percentage of ambiguous mappings, overall alignment rate.
(in jobs/Job_X/STAT/bowtie2.log)

5.2.3 Identification/Annotation:

- BLAT: As specified in the header, BLAT reports various alignment statistics (percentage of identity, alignment length, number of mismatches, number of gap opening, E value and Bit score)
(in jobs/Job_X/IDENTIFICATION/blat_transcripts/DB/blat_DB_summary.txt, "DB" being the name of the database you specified).
For more information see:
<http://genome.ucsc.edu/goldenpath/help/blatSpec.html>
- HMMER: As specified in the header, HMMER reports various profile match statistics (such as E value and score)
(in jobs/Job_X/IDENTIFICATION/hmmer_transcripts/DB/hmmer_results.tab.txt, "DB" being the name of profile database you specified).
For more information see:
<http://hmmer.janelia.org/>
- BLAST+: All common BLAST statistics (such as percentage of identity, E value and Bit Score) are available in an xml formatted file.
(in jobs/Job_X/IDENTIFICATION/blastplus_par/blast_out/total.xml).

Finally a file summarizing the BLAT and BLAST results is available in jobs/Job_X/STAT/total_annotation.csv and can be open in excel.

5.2.4 Expression

- RSEM: As specified in the headers, effective length of the transcripts, expected counts, TPM, RPKM are provided both at the isoform and gene level.
(in jobs/Job_X/EXPRESSION/RSEM/RSEM.isoforms.results and RSEM.genes.results).
For more information see:
<http://deweylab.biostat.wisc.edu/rsem/>

- eXpress: As specified in the header, expression statistics such as total counts, effective counts, TPM and RPKM are provided for each isoforms.
(in jobs/Job_X/EXPRESSION/express/results.xprs).
For more information see:
<http://bio.math.berkeley.edu/eXpress/manual.html>

6 Other useful programs:

Part of Truفا's output can be visualized with programs such as :

- Blast2GO
<http://www.blast2go.com/b2ghome>
- IGV
<http://www.broadinstitute.org/igv/>
- Tablet
<http://bioinf.scri.ac.uk/tablet/>
- RNAseq Viewer
<http://bioinfo.au.tsinghua.edu.cn/software/RNAseqViewer/>

Other programs for RNA-seq analysis can be found at:

- <http://omictools.com/>
- [Wikipedia: RNA-seq tools](#)

7 Quick Start:

7.1 Launch a Job

This is a quick demonstration on how the web server can be used.

7.1.1 Upload demo files



- In the Upload area, **Browse** your computer for the reads file “reads_left.fq.tar.gz”
- Before clicking **Send**, select **compressed reads file** in the dropdown menu (instead of **undefined**)
- Once successfully uploaded, do the same for the reads file “reads_right.fq.tar.gz”
- Once both files have been uploaded, you are ready to “start a job”.

Specifying the file type in the dropdown menu is mandatory and will help you by filtering the correct inputs for each steps of the analysis.

7.1.2 Start the job



- Go to **Start a Job**
- Select **Paired-end reads (2 fastq files)**
- Specify the 2 input files: “reads_left.fq.tar.gz” for **left reads** and “reads_right.fq.tar.gz” for **right reads**
- Go down to **RNA-seq steps** section and to the **Reads cleaning** tab.
- Check the boxes **Duplicated reads** and **Quality trimming** in the Prinseq section
- Check the box **FastQC** in the **Post-cleaning quality control** section
- Then go to the **Reads Assembly and Mapping** tab
- Check the boxes **Assemble with Trinity** and **Align reads against contigs with Bowtie2**
- Then go to the **Contigs identification** tab
- Go down to **Blast2GO searches** section
- Check the box **Blast2GO**
- Then go to the **Expression quantification** tab
- Check **eXpress**
- Then go down to **Launching Analysis** and hit the **Start** button.

A message confirming that the Job has been sent should appear. You can go back to **Home** in the main menu and here click on the corresponding Job in the Job list. You should be able to see the current state of your job there.

7.2 Check output

Once your job has been completed, here are some ways to check the generated outputs.



- Go to **File Manager** in the main menu
- Put your username and password
- Browse the **Jobs** folder for the corresponding Job name

7.2.1 Check reads quality:



- Go to **STAT** and then **prinseq_report**

Here you should have 4 html files (2 for each reads file, one report after the removal of the duplicates and one report after the trimming process) that you can download and then open locally with your internet browser. Information about Prinseq statistics can be found in the “quality control” section of Prinseq manual:

<http://prinseq.sourceforge.net/manual.html#QC>



- Go to **STAT** and then **fastqc_report**

Here you should have 2 zip files (1 for each reads file). Download them, extract them and open locally the html files with your internet browser. Information about the FastQC output can be found here:

<http://www.bioinformatics.babraham.ac.uk/projects/fastqc/>

7.2.2 Check the assembly:



- Go to **ASSEMBLY MAPPING** and then **trinity**

Here you will find the assembly file **Trinity.fasta**



- Go to **STAT** and then **assembly_qc**

Here you will find text files displaying statistics related to your assembly which will help you to determine the quality of your assembly.

7.2.3 Check reads mapping against transcripts



- Go to **ASSEMBLY MAPPING** and then **bowtie2**
- Download the file **aligned_reads.bam** and **aligned_reads.bam.bai**

You can use these two files to visualize the alignments using programs such as Tablet <http://bioinf.scri.ac.uk/tablet/>.

7.2.4 Check the Blast2GO results:



- Go to **IDENTIFICATION** and then **b2go**
- Download the file **out_b2go.dat**
- Open it by starting Blast2GO program locally and going to **File, Load B2GO-Project (.dat)**

Blast2GO can be obtained from <http://www.blast2go.com/b2ghome>

7.2.5 Check the eXpress results:



- Go to **EXPRESSION** and then **express**

Here you will find two text files (“params.xprs” and “results.xprs”). The “results.xprs” file is giving the expression quantification for all isoforms produced by Trinity. More information about eXpress outputs can be found at <http://bio.math.berkeley.edu/eXpress/manual.html>

8 Running an analysis:

This explains in more details each parts which can be found in the “Start a Job” page of the web server. This manual is far from pretending to replace the manual of each program which are implemented in TRUFA. Consequently, the manuals of each of these programs should be studied by the user in order to tune them and understand their output properly.

8.1 Type of input:

Essentially, all the steps of the RNA-seq pipeline on the web server can be realized with only reads files (1 or 2) as input. For the purpose of testing and parameters tuning, later steps can be as well directly performed with an assembly (fasta file with already assembled reads).

First, you should specify and upload the necessary input files:

- For Reads cleaning: You will need one or two fastq files (i.e. single or paired-end reads files)
- For Assembly: You will need one or two fastq files (i.e. single or paired-end reads files)
- For Identification, you will need either: a fasta file with the reads already assembled OR 1 or 2 fastq reads files and generate an assembly.
- For Reads mapping: 1 fasta file with the reads already assembled and the corresponding reads in 1 or 2 fastq files OR 1 or 2 fastq reads files and input or generate an assembly.
- For Expression quantification: 1 fasta file with the reads already assembled and the corresponding reads in 1 or 2 fastq files OR 1 or 2 fastq reads files and input or generate an assembly (see How to generate an assembly).

Then, by specifying in the **Run a Job** part the type of input and the input files, different possibilities will be available to you, depending on the type and number of input files you specified.

8.2 Cleaning step:

In this well, you will find programs to perform reads cleaning prior to the assembly.

- **FastQC:** Will give an idea on how to clean your reads. Are your reads of poor quality ? Got hexamer priming bias ? Got contaminants like adapters sequences in my reads ?

[More info on FASTQC homepage](#)

- **Cutadapt:** Allows you to remove the adapter sequences you might have found during the FASTQC quality control. Adapter sequences must be specified further down in the “**Cutadapt options**”.

[More on Cutadapt](#)

- **Prinseq:** Performs duplicate removal and reads trimming. Lot's of options can be tuned further down the “Start a Job” page in “**Duplication options**” and “**Trimming options**”. At this stage, removing duplicates is mainly to gain computation time during the assembly process. For expression quantification purposes, it seems more recommended to mark duplicates using programs such as Picard.

[More info on Prinseq homepage](#)

- **Blat:** performs homology searches but much faster than Blast and can be therefore used to search homologies between a large set of reads and a database, in order to remove potential contaminants before assembly. Custom databases of potential contaminants can be added by **Uploading** them as “**nucleotide seqs database for blasting**” in fasta format. Then the corresponding database can be checked in the “**Nucleotide db**” button in the “**Blast against potential contaminants**” area. Blat parameters cannot be tuned on TRUFA yet and is used with default parameters. Custom scripts are automatically removing for the next steps of the analysis the reads with at least one blat hit for any of the selected database

[More info on Blat](#)

8.3 Assembly and Mapping step:

In this well, the programs will perform assembly and/or mapping of the reads. If Cutadapt, Prinseq or Blat actions have been checked in the “Cleaning

step” well, then the output of this cleaning step will be used as input for the assembly and mapping steps. If no cleaning options have been selected, raw reads will be used for the assembly/mapping.

- **Trinity:** performs de novo assembly of the reads. Several options are important to tune, such as the minimum contig length and the strand specificity if required. This can be tuned in the “**Trinity options**”.
[More info on Trinity](#)
- **Assembly quality checks:** provide various statistics helping to judge the quality of the assembly. This includes the [CEGMA](#) program, [blast searches](#) and [Trinity scripts](#) and custom scripts for statistics and plots.
- **Bowtie2:** align the reads back to the transcripts. on top of being necessary for the expression quantification, this mapping and the amount of reads which align back are an indicator of the quality of the assembly. It’s important to verify that the insert size of your library is comprise between the minimum and the maximum insert size in “**Bowtie2 options**”.
[More info on Bowtie2](#)

8.4 Identification step:

This part of the pipeline aims at annotating the transcripts newly assembled and link nucleotide sequences to a particular biological function.

- **Blat:** Again, Blat is here to perform fast homology searches against classical databases such as NCBI nr or Uniref90. Blat parameters can not be tuned yet and are set to default. Custom databases can be added by **Uploading** them as “**nucleotide (or protein) seqs database for blasting**” in fasta format. Then the corresponding database can be checked in the respective button (“nucleotide” or “protein” db) in the “**Custom Blat searches**” area.
[More info on Blat](#)
- **HMMER:** performs as well homology searches but is able to detect more remote homologs with HMM protein profiles. The HMMER search can be performed on the whole PFAMA database or against one/multiple custom hmm profiles previously generated (TRUFA does

not generate HMM profile files) and **Uploaded** as “**protein hmm profile for HMMER**” in .hmm format.

[More info on HMMER](#)

- **Blast+ against nr and B2GO:** To get the GO terms annotations for each transcripts, B2GO first needs the results of a blast search against the NCBI nr database. The Blast+ job is the most computation intensive from TRUFA and can in certain cases stop before the end of the search is reached. If you notice that B2GO is not producing any output, please contact us, this could be the issue. [More info on Blast+](#)
[More info on Blast2GO](#)

8.5 Expression quantification step:

This step allows the user to produce classical expression quantification outputs with values such as TPM, FPKM/RPKM or counts values useful to further test for differential expression (Differential expression analysis is not yet available on TRUFA).

9 Security and privacy note:

TRUFA is using the official certificate from the National Research and Academic Network (RedIRIS, <http://www.rediris.es/>) to approve for the security of its https connection. Data uploaded to the webserver are only accessible to the user who uploaded it and to the developers Etienne Kornobis, Luis Cabellos and Fernando Aguilar only for support purposes.

References

Jeffrey A Martin and Zhong Wang. Next-generation transcriptome assembly. *Nature reviews. Genetics*, 12(10):671–82, October 2011. ISSN 1471-0064. doi: 10.1038/nrg3068. URL <http://www.ncbi.nlm.nih.gov/pubmed/21897427>.

Jochen B W Wolf. Principles of transcriptome analysis and gene expression quantification: an RNA-seq tutorial. *Molecular ecology resources*, 13:559–572, April 2013. ISSN 1755-0998. doi: 10.1111/1755-0998.12109. URL <http://www.ncbi.nlm.nih.gov/pubmed/23621713>.